

Protective Overall CK-1

Categorization

Protection class: Type 3.

Product description

The CK1 consists of a body piece, long sleeves, long trouser legs and a hood, attached with an overlock-seam and a zipper cover strip. The zipper covers the chin and the neck areas. The cut of the coverall is designed for more mobility. The coverall is equipped with an overlock-seam bordered elastic strap at the face, arm, leg, and waist. Material: SMMS.

Material

SMMS 60g

Applied standards

EN 14126:2003, EN 14325:2001,
ISO/FDIS 16603 and ISO/FDIS 16604,
ISO/DIS 22611 and ISO/DIS 22622,
EN 340,
Type 3:EN466

Application

This coverall offers protection against dangerous substances and contamination. It can be used, depending on the circumstances and degree of the toxicity, against air carried particles and nontoxic splashes and spray



For general isolation of outpatient/ward/laboratory in medical institutions



For civil protection such as security inspection/disinfection/sentry box