

Isolation Gown

Isolation Gown is made of non woven fabric and SMS fabric, soft and skin-friendly.

Used in hospital, laboratory to prevent and isolate dust particle, alcohol, blood, bacteria and virus invading.

Material: non woven fabric (classic gown), MS fabric (premium gown), reinforced gown.

Cuff: knitted cuff, elastic cuff.

With 2 waist belts, 2 neck belts (or velcro belts)

Color: white, green, blue, pink, yellow. etc.

Size: 115x137cm, 120x140cm, 130x150cm, 140x60cm, S, M, L, XL, XXL

Packing: 1 piece in a sterile pouch, 50 pcs/carton

Sterile by EO gas if need.

